

BS EN 16091:2011
BS 2000-595:2011



BSI Standards Publication

**Liquid petroleum products —
Middle distillates and fatty
acid methyl ester (FAME) fuels
and blends — Determination
of oxidation stability by rapid
small scale oxidation method**

bsi.

...making excellence a habit.™

EUROPEAN STANDARD

EN 16091

NORME EUROPÉENNE

EUROPÄISCHE NORM

November 2011

ICS 75.160.20

English Version

Liquid petroleum products - Middle distillates and fatty acid
methyl ester (FAME) fuels and blends - Determination of
oxidation stability by rapid small scale oxidation method

Produits pétroliers liquides - Distillats moyens, esters
méthyliques d'acides gras (EMAG) et leurs mélanges -
Détermination de la stabilité à l'oxydation par méthode
d'oxydation accélérée à petite échelle

Flüssige Mineralölerzeugnisse - Mitteldestillat- und
Fettsäuremethylesterkraftstoffe und Mischungen -
Bestimmung der Oxidationsstabilität mit beschleunigtem
Verfahren und kleiner Probenmenge

This European Standard was approved by CEN on 15 October 2011.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN-CENELEC Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

Management Centre: Avenue Marnix 17, B-1000 Brussels

Contents		Page
Foreword.....		3
1	Scope	4
2	Normative references	4
3	Terms and definitions	4
4	Principle	4
5	Reagents and materials	4
6	Apparatus	5
6.1	Automatically controlled oxidation tester, comprising an oxidation reaction vessel containing:	5
6.2	Pipette, capable of delivering 5,0 ml ± 0,1 ml.	5
6.3	O-ring seals, see A.6.	5
7	Sampling and sample handling.....	5
8	Performance check of the apparatus.....	5
9	Preparation of the apparatus.....	6
10	Procedure	6
11	Expression of results	7
12	Precision	7
12.1	General.....	7
12.2	Repeatability, <i>r</i>	7
12.3	Reproducibility, <i>R</i>	7
13	Test report	7
Annex A (normative) Apparatus description.....		8
A.1	General requirements.....	8
A.2	Assembly of apparatus	8
A.3	Reaction vessel.....	9
A.4	Screw cap	9
A.5	Electric heater	9
A.6	O-ring seals	9
A.7	Valves.....	9
A.8	Pressure sensor.....	9
A.9	Temperature sensor	9
A.10	Connecting pipes.....	9
A.11	Cooling fan	10
Annex B (informative) Calibration process.....		11
Annex C (informative) Determination of induction period		12
Bibliography		13

Foreword

This document (EN 16091:2011) has been prepared by Technical Committee CEN/TC 19 “Gaseous and Liquid fuels, lubricants and related products of petroleum, synthetic and biological origin”, the secretariat of which is held by NEN.

This European Standard shall be given the status of a national standard, either by publication of an identical text or by endorsement, at the latest by May 2012, and conflicting national standards shall be withdrawn at the latest by May 2012.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN [and/or CENELEC] shall not be held responsible for identifying any or all such patent rights.

According to the CEN/CENELEC Internal Regulations, the national standards organizations of the following countries are bound to implement this European Standard: Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and the United Kingdom.

1 Scope

This European Standard specifies a method for the determination of the oxidation stability of middle distillate fuels, fatty acid methyl ester (FAME) fuel and blends thereof, under accelerated conditions, by measuring the induction period to the specified breakpoint in a reaction vessel charged with the sample and oxygen.

NOTE 1 For the purposes of this European Standard, the term “% (V/V)” is used to represent the volume fraction (φ).

NOTE 2 The induction period is used as an indication for the resistance of middle distillates, fatty acid methyl ester (FAME) fuels and blends thereof against oxidation. It should be recognized, however, that this correlation can vary markedly under different conditions with different FAMEs and diesel fuel blends.

NOTE 3 The presence of ignition improvers may lead to lower oxidation stability results determined by this method. It has for instance been observed that the addition of 2-ethyl hexyl nitrate (2EHN) can reduce the measured oxidation stability values.

2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN ISO 3170, *Petroleum liquids — Manual sampling (ISO 3170:2004)*

EN ISO 3171, *Petroleum liquids — Automatic pipeline sampling (ISO 3171:1988)*

3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

3.1 breakpoint

point in the pressure-time curve, where the pressure in the vessel has dropped 10 % below the maximum pressure achieved in the current test run

3.2 induction period

time elapsed between starting the heating procedure of the reaction vessel, charged with sample and oxygen, and the breakpoint at 140 °C

4 Principle

At ambient temperature, a known volume of a sample is placed in a reaction vessel charged with oxygen to a pressure of 700 kPa \pm 5 kPa. The reaction vessel is heated to 140 °C. The pressure in the vessel drops as the oxygen is consumed during the oxidation of the sample. The pressure in the vessel is recorded at intervals of 1s until the breakpoint is reached. The elapsed time from start to the breakpoint is the induction period at the test temperature of 140 °C \pm 0,5 °C.

5 Reagents and materials

5.1 Cleaning solvent, for the removal of oxidation residues from the test vessel, of suitable purity to leave no residue on the apparatus.

NOTE Commercially available ethanol of approx. 95 % (V/V) purity was found to be suitable.

5.2 Oxygen, extra dry (< 5 mg/kg water), commercially available, with a purity of not less than 99,6 %.

5.3 Cleaning tissues, lint-free for cleaning sensitive galvanic coated surfaces without scratching.

5.4 Verification fluid, hydrocarbon fuel (B0).

NOTE Any hydrocarbon fuel with sufficient stability and a known induction period may be used as verification fluid. In general, verification fluid with a certified induction period is available from the manufacturer of the apparatus.

6 Apparatus

6.1 Automatically controlled oxidation tester, comprising an oxidation reaction vessel containing:

- test sample cup capable of being rapidly heated;
- pressure sensor capable of measuring pressures of 1 kPa up to 2 000 kPa with a resolution of 1 kPa;
- temperature sensor with a resolution of 0,1 °C;
- pressure and temperature recording of 1 second.

The oxidation reaction vessel shall be fitted with filling and relief valves, meant to release the pressure, and a fan to cool the vessel from test to ambient temperature.

The requirements for the apparatus are described in detail in Annex A.

6.2 Pipette, capable of delivering 5,0 ml ± 0,1 ml.

6.3 O-ring seals, see A.6.

7 Sampling and sample handling

Unless otherwise specified, samples shall be taken as described in EN ISO 3170 or EN ISO 3171 and/or in accordance with the requirements of national standards or regulations for the sampling of the product under test.

Collect and store samples in an opaque container to minimize exposure to light.

8 Performance check of the apparatus

8.1 Verify the correct operation of the equipment by following the procedure in section 10 using the verification fluid (5.4).

8.2 If the induction period of the verification fluid does not meet specified time, repeat 8.1. If it fails a second time, refer to the manufacturer's instructions.

8.3 Pressure and temperature are critical parameters. Proper calibration of the respective sensors (see 6.1) is therefore important. Recalibrate the tester at least every 12 months for correct temperature and pressure readings.

NOTE In view of general good laboratory practice and harmonization of test methods, calibration as in Annex B is strongly recommended. As the calibration procedure for a pressure sensor is complicated and difficult to be executed by an untrained lab technician, one may (in most cases sufficiently) rely on the verification fluid.

9 Preparation of the apparatus

9.1 Remove the previous test portion by means of a pipette or similar device.

9.2 Remove the used "O-ring" seal and discard.

NOTE To avoid contamination of the new test, it is necessary to discard the used "O-ring" seal, because it may be soaked with oxidation products from the previous test.

9.3 Wipe the sample cup, the seal groove and the cover with lint-free cleaning tissue (5.3) soaked with cleaning solvent (5.1) until free of gum or other oxidation residues.

9.4 Allow the test sample cup and cover to dry in air and visually inspect for cleanliness.

NOTE Compressed air is generally unsuitable to accelerate the evaporation of solvent because it can contain traces of oil that may influence the next test.

9.5 Insert a new "O-ring" seal (6.3).

10 Procedure

10.1 Bring the reaction vessel and the sample to be tested to room temperature.

10.2 Using a pipette (6.2), place $5\text{ ml} \pm 0,1\text{ml}$ of the sample into the reaction vessel sample cup.

10.3 Cover the sample-cup with the screw cap (see A.2) and close the reaction vessel.

10.4 Pre-flush the reaction vessel with oxygen. Introduce oxygen (5.2) into the vessel until a pressure of $700\text{ kPa} \pm 5\text{ kPa}$ is attained and stabilized over 20 sec.

10.5 Start the heater with no delay between charging with oxygen and starting the test.

10.6 The test temperature of $140\text{ }^{\circ}\text{C} \pm 0,5\text{ }^{\circ}\text{C}$ shall be reached after $270\text{ s} \pm 30\text{ s}$.

10.7 If, during the initial 5 min of the test, a steady drop in pressure is observed, discontinue the test and discard the test specimen.

The leakage rate at $140\text{ }^{\circ}\text{C}$ with empty sample cup shall not exceed a value of 2 kPa/h . If the leakage rate shows an increase, check the following components:

- O-ring for damage or residues of samples,
- surface of sample cup for damage, and
- sample cup for sample residues.

Contact the manufacturer to resolve leakage problems from other parts of the instrument.

10.8 The apparatus automatically records the temperature (to the nearest $0,1\text{ }^{\circ}\text{C}$) and pressure (to the nearest 1 kPa) of the oxidation vessel continuously.

10.9 The apparatus terminates the test when the pressure readings show a 10 % drop from the maximum observed pressure (breakpoint) and records the time elapsed from the start of the test to the breakpoint (induction period). An example is given in Annex C.

10.10 The apparatus automatically records the average temperature to the nearest $0,1\text{ }^{\circ}\text{C}$ as the temperature of the test.

NOTE The apparatus will automatically switch on the fan to cool the reaction vessel to approximately room temperature. When the apparatus has cooled, the pressure is slowly released, at a rate not exceeding 345 kPa/min.

10.11 When the pressure release has finished, open the apparatus and clean according to Clause 9.

11 Expression of results

Report the induction period expressed in min., rounded to the nearest 0,01 min.

12 Precision

12.1 General

The precision statements were developed in Round Robin tests performed in 2008 and 2009, in accordance with EN ISO 4259 [1]. This method was tested for FAME, B0, B5, B7, B10, and B30 samples, the induction period measured being between 22 min and 215 min.

12.2 Repeatability, r

The difference between two test results, obtained by the same operator with the same apparatus under constant operating conditions on identical test material, would in the long run, in the normal and correct operation of the test method, exceed the values calculated from the following equation only in one case in twenty:

$$r = 0,028\,8X + 0,496\,5 \quad (1)$$

where

X represents the mean of the two results expressed in min, rounded to the nearest 0,01 min.

12.3 Reproducibility, R

The difference between two single and independent test results, obtained by different operators working in different laboratories on identical test material, would in the long term, in the normal and correct operation of the test method, exceed the values calculated from the following equation only in one case in twenty:

$$R = 0,086\,3X + 1,377\,2 \quad (2)$$

where

X represents the mean of the two results expressed in min, rounded to the nearest 0,01 min.

13 Test report

The test report shall contain at least the following information:

- a) a reference to this European Standard;
- b) the type and complete identification of the product tested;
- c) the result of the test (see 11);
- d) any deviation, by agreement or otherwise, from the procedure specified;
- e) the date of the test.

Annex A
(normative)


Apparatus description

A.1 General requirements

All equipment in contact with the sample and oxygen shall be made of corrosion-resistant materials and the reaction vessel shall have an over-temperature safeguard.

A.2 Assembly of apparatus

The assembly shall be according to Figure A.1.



Key	
1	screw cap (A.4)
2	O-ring seal(A.6)
3	reaction vessel (A.3)
4	oxygen supply (5.2)
5	temperature sensor (A.9)
6	electric heater (A.5)
7	cooling (A.11)

Figure A.1 — Reaction vessel for accelerated oxidation

A.3 Reaction vessel

Inner diameter 47 mm and internal volume without sample of 20 ml. Made of aluminium alloy or other corrosion-resistant metal of similar thermal conductivity, with a cylindrical depression (test cup), and with a hole in the side to accommodate a temperature sensor (A.9).

The interior surface of the test cup is covered with galvanic coating and smoothly polished.

A.4 Screw cap

Made of aluminium alloy or other corrosion-resistant metal of similar thermal conductivity. Interior surface of the screw cap covered with galvanic coating, smoothly polished.

A.5 Electric heater

Power rating 500 W; attached to the bottom of the test cup in a manner that provides efficient transfer of heat. The heater control shall be capable of maintaining the test cup temperature within $\pm 0,5$ °C for test temperatures up to and including 140 °C.

A.6 O-ring seals

Made of material being resistant to oxygen, middle distillate fuel components, FAME and heat, typically fluoro-elastomer FPM/FKM (commonly known as Viton®¹), coated with PTFE.

A.7 Valves

Solenoid with small dead volume and orifice sizes and short cycling rates.

A.8 Pressure sensor

Absolute or relative pressure transducer with a pressure range of 0 kPa to 2 000 kPa and a typical sensitivity of between (1 to 2,5) mV/kPa. Capable of being operated at temperatures up to 150 °C and calibrated and/or verified against an instrument certified by an authorized certification body.

A.9 Temperature sensor

Platinum thermometer capable of reading up to 200 °C $\pm 0,1$ °C, calibrated and/or verified against an instrument certified by an authorized certification body.

A.10 Connecting pipes

Typical inner diameter: 0,5 mm to 1 mm.

¹ Viton is an example of a suitable product available commercially. This information is given for the convenience of users of this European Standard and does not constitute an endorsement by CEN of this product.

A.11 Cooling fan

Capable of cooling, the reaction vessel from the test temperature to ambient by applying a stream of air to the outside of the reaction vessel.

Annex B (informative)

Calibration process

B.1 Calibration of temperature indicator.

Calibrate the temperature sensor to the nearest 0,1 °C according to the manufacturer's instructions.

B.2 Calibration of pressure detector.


Calibrate the pressure sensor to the nearest 10 kPa according to the manufacturer's instructions.

B.3 Verify, at least every 3 months, that the heater is operating properly and the reaction vessel is reaching $140\text{ °C} \pm 0,5\text{ °C}$ within $270\text{ s} \pm 30\text{ s}$.

Annex C
(informative)

Determination of induction period

Figure C.1 presents an example of how the breakpoint and induction period of a sample are determined. The point B indicates the 10% drop from maximum pressure, $0,9 P_{\text{max}}$, which is the breakpoint and defines the induction period (C).



- Key**

A maximum observed pressure, P_{max}

B breakpoint

C induction period (min)
- D test time (min)

E observed pressure (kPa)

Figure C.1 — Typical pressure/time curve

Bibliography

- [1] EN ISO 4259, *Petroleum products — Determination and application of precision data in relation to methods of test (ISO 4259:2006)*

Buying Parts of BS 2000

Orders for BS 2000 publications should be addressed to either:

Energy Institute – Library and Information Service

61 New Cavendish Street
London
W1G 7AR

Tel: +44 (0)20 7467 7100
Fax: +44 (0)20 7255 1472
www.energyinst.org.uk

Order standards securely via:
www.energyinstpubs.org.uk

or:

British Standards Institution – Customer Services

389 Chiswick High Road
London
W4 4AL

Tel: +44 (0)20 8996 9001
Fax: +44 (0)20 8996 7001
www.bsigroup.com

Order hard copy standards securely via:
www.bsigroup.com/shop

Copyright

Copyright exists in all BS 2000 publications. No part of this publication may be reproduced in any form without the prior permission in writing of BSI and the Energy Institute. Enquiries about copyright should be made to the Secretary of PTI/13 at the Energy Institute.



British Standards Institution (BSI)

BSI is the national body responsible for preparing British Standards and other standards-related publications, information and services.

BSI is incorporated by Royal Charter. British Standards and other standardization products are published by BSI Standards Limited.

About us

We bring together business, industry, government, consumers, innovators and others to shape their combined experience and expertise into standards-based solutions.

The knowledge embodied in our standards has been carefully assembled in a dependable format and refined through our open consultation process. Organizations of all sizes and across all sectors choose standards to help them achieve their goals.

Information on standards

We can provide you with the knowledge that your organization needs to succeed. Find out more about British Standards by visiting our website at bsigroup.com/standards or contacting our Customer Services team or Knowledge Centre.

Buying standards

You can buy and download PDF versions of BSI publications, including British and adopted European and international standards, through our website at bsigroup.com/shop, where hard copies can also be purchased.

If you need international and foreign standards from other Standards Development Organizations, hard copies can be ordered from our Customer Services team.

Subscriptions

Our range of subscription services are designed to make using standards easier for you. For further information on our subscription products go to bsigroup.com/subscriptions.

With **British Standards Online (BSOL)** you'll have instant access to over 55,000 British and adopted European and international standards from your desktop. It's available 24/7 and is refreshed daily so you'll always be up to date.

You can keep in touch with standards developments and receive substantial discounts on the purchase price of standards, both in single copy and subscription format, by becoming a **BSI Subscribing Member**.

PLUS is an updating service exclusive to BSI Subscribing Members. You will automatically receive the latest hard copy of your standards when they're revised or replaced.

To find out more about becoming a BSI Subscribing Member and the benefits of membership, please visit bsigroup.com/shop.

With a **Multi-User Network Licence (MUNL)** you are able to host standards publications on your intranet. Licences can cover as few or as many users as you wish. With updates supplied as soon as they're available, you can be sure your documentation is current. For further information, email bsmusales@bsigroup.com.

BSI Group Headquarters

389 Chiswick High Road London W4 4AL UK

Revisions

Our British Standards and other publications are updated by amendment or revision.

We continually improve the quality of our products and services to benefit your business. If you find an inaccuracy or ambiguity within a British Standard or other BSI publication please inform the Knowledge Centre.

Copyright

All the data, software and documentation set out in all British Standards and other BSI publications are the property of and copyrighted by BSI, or some person or entity that owns copyright in the information used (such as the international standardization bodies) and has formally licensed such information to BSI for commercial publication and use. Except as permitted under the Copyright, Designs and Patents Act 1988 no extract may be reproduced, stored in a retrieval system or transmitted in any form or by any means – electronic, photocopying, recording or otherwise – without prior written permission from BSI. Details and advice can be obtained from the Copyright & Licensing Department.

Useful Contacts:

Customer Services

Tel: +44 845 086 9001

Email (orders): orders@bsigroup.com

Email (enquiries): cservices@bsigroup.com

Subscriptions

Tel: +44 845 086 9001

Email: subscriptions@bsigroup.com

Knowledge Centre

Tel: +44 20 8996 7004

Email: knowledgecentre@bsigroup.com

Copyright & Licensing

Tel: +44 20 8996 7070

Email: copyright@bsigroup.com



...making excellence a habit.™