



BSI Standards Publication

Plastics — Determination of caprolactam and its cyclic and linear oligomers by HPLC

National foreword

This British Standard is the UK implementation of ISO 15033:2018. It supersedes BS ISO 15033:2007, which is withdrawn.

The UK participation in its preparation was entrusted to Technical Committee PRI/21, Testing of plastics.

A list of organizations represented on this committee can be obtained on request to its secretary.

This publication does not purport to include all the necessary provisions of a contract. Users are responsible for its correct application.

© The British Standards Institution 2018
Published by BSI Standards Limited 2018

ISBN 978 0 580 94256 3

ICS 83.080.01

Compliance with a British Standard cannot confer immunity from legal obligations.

This British Standard was published under the authority of the Standards Policy and Strategy Committee on 30 June 2018.

Amendments/corrigenda issued since publication

Date	Text affected
------	---------------

**INTERNATIONAL
STANDARD**

**ISO
15033**

Third edition
2018-06-15

**Plastics — Determination of
caprolactam and its cyclic and linear
oligomers by HPLC**

*Plastiques — Détermination du caprolactame et de ses oligomères
cycliques et linéaires par CLHP*



Reference number
ISO 15033:2018(E)

© ISO 2018



COPYRIGHT PROTECTED DOCUMENT

© ISO 2018, Published in Switzerland

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
Ch. de Blandonnet 8 • CP 401
CH-1214 Vernier, Geneva, Switzerland
Tel. +41 22 749 01 11
Fax +41 22 749 09 47
copyright@iso.org
www.iso.org

Contents		Page
Foreword		iv
Introduction		v
1 Scope		1
2 Normative references		1
3 Terms and definitions		1
4 Principle		1
5 Reagents		1
6 Apparatus		3
7 Test sample		4
8 Procedure		4
8.1 Calibration		4
8.1.1 Method A.....		4
8.1.2 Method B.....		5
8.2 Determination		5
8.2.1 Method A.....		5
8.2.2 Method B.....		7
9 Calculations		7
9.1 Calculation of calibration factors		7
9.1.1 Calibration factors for cyclic oligomers		7
9.1.2 Calibration factors for linear oligomers.....		8
9.2 Test sample.....		8
9.2.1 Calculation of cyclic oligomer concentrations.....		8
9.2.2 Calculation of linear oligomer concentrations.....		9
10 Precision		9
11 Test report		10
Annex A (normative) HPLC parameters and injector programme — Using eluent A		11
Annex B (normative) HPLC parameters and injector programme — Using eluent B		12
Annex C (informative) Schematic diagram of HPLC apparatus		13
Bibliography		14

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

The procedures used to develop this document and those intended for its further maintenance are described in the ISO/IEC Directives, Part 1. In particular the different approval criteria needed for the different types of ISO documents should be noted. This document was drafted in accordance with the editorial rules of the ISO/IEC Directives, Part 2 (see www.iso.org/directives).

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights. Details of any patent rights identified during the development of the document will be in the Introduction and/or on the ISO list of patent declarations received (see www.iso.org/patents).

Any trade name used in this document is information given for the convenience of users and does not constitute an endorsement.

For an explanation on the voluntary nature of standards, the meaning of ISO specific terms and expressions related to conformity assessment, as well as information about ISO's adherence to the World Trade Organization (WTO) principles in the Technical Barriers to Trade (TBT) see the following URL: www.iso.org/iso/foreword.html.

This document was prepared by Technical Committee ISO/TC 61, *Plastics*, Subcommittee SC 5, *Physical-chemical properties*.

This third edition cancels and replaces the second edition ([ISO 15033:2007](http://www.iso.org/iso/15033:2007)), which has been technically revised to add method B.

Introduction

The method A specified in this document can be used for HPLC determination of the cyclic oligomers of caprolactam ($n = 1$) up to and including the hexamer ($n = 6$), using UV detection. If desired, after post-column reaction of the primary amine with 1,2-phthalic dicarboxaldehyde, online determination of the linear oligomers up to and including the hexamer can also be carried out.

The determination is not quantitative for oligomers higher than the hexamer ($n > 6$). In the determination of cyclic oligomers the sensitivity for the tetramer and higher oligomers is constant, which means that calibration should take place up to and including the tetramer ($n = 4$).

The linear oligomers are determined by the fluorescence of the iso-indole group, which is a product of the reaction between the primary amino group, 1,2-phthalic dicarboxaldehyde and 3 mercaptopropionic acid. The calibration with the linear oligomers should be carried out up to and including the hexamer ($n = 6$).

The method B included in this document is intended for the determination of caprolactam and its cyclic dimer only, following the same principle and using the same equipment as method A, but significantly faster.

Plastics — Determination of caprolactam and its cyclic and linear oligomers by HPLC

SAFETY STATEMENT — Persons using this document should be familiar with normal laboratory practice, if applicable. This document does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user to establish appropriate safety and health practices and to determine the applicability of any other restrictions.

1 Scope

This document describes a high-performance liquid chromatography (HPLC) method for determining the concentrations of cyclic oligomers of caprolactam, from 0,01 % by mass upwards, and linear oligomers of caprolactam, from 5 mg/kg upwards, both up to and including the hexamer of caprolactam ($n = 6$), in samples of polyamide 6, caprolactam and mixtures of rearrangement products in water.

A second, significantly faster, HPLC method is included for determination of caprolactam and its cyclic dimer, based on the same principle and using the same equipment as the first method.

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

[ISO 472](#), *Plastics — Vocabulary*

3 Terms and definitions

For the purposes of this document, the terms and definitions given in [ISO 472](#) apply.

ISO and IEC maintain terminological databases for use in standardization at the following addresses:

- IEC Electropedia: available at <http://www.electropedia.org/>
- ISO Online browsing platform: available at <https://www.iso.org/obp>

4 Principle

A test sample is dissolved in, or diluted with, formic acid and the oligomers separated in the presence of a low-pH mobile phase using a column filled with reversed-phase packing material. The cyclic oligomers are detected by UV absorption at 200 nm. If desired, the linear oligomers can be detected by fluorescence after post-column reaction of the primary amino group with 1,2-phthalic dicarboxaldehyde and 3-mercaptopropionic acid. The concentrations are calculated by comparison of the measured values with those of calibration solutions.

5 Reagents

During the analysis, unless otherwise stated, use only reagents of recognized analytical grade.

5.1 Water, ultrapure or double-distilled.

5.2 Phosphoric acid, 85 % by mass.

5.23 Post-column derivatization reagent.

Dissolve 76 g of sodium tetraborate decahydrate (5.18) and 6 g of sodium hydroxide (5.19) in 2 l of water (5.1). Dissolve 1,6 g of 1,2-phthalic dicarboxaldehyde (5.20) in 40 ml of methanol (5.21) and add this solution to the sodium tetraborate decahydrate reagent. Add 1,5 ml of 3-mercaptopropionic acid (5.22) and mix well.

The stability of the post-column derivatization reagent is limited. Do not keep for longer than 3 days.

NOTE The cyclic dimer, the cyclic trimer and the mixture of cyclic oligomers of caprolactam can be isolated from a methanol extract of polyamide 6 (PA6) by preparative HPLC, using the HPLC methods described here. The purity of the dimer and the possible presence of other oligomers can be checked using the methods described in this document.

6 Apparatus

6.1 HPLC equipment, having the following specifications:

- **Eluent pump**, including mixer, damper and manometric module, giving an eluent flow rate of 0,51 ml/min and a pressure drop of approximately 10 MPa.
- **Injector**, e.g. an auto-sampler capable of 1 µl to 250 µl injections, equipped to carry out a variable injector programme (see Annex A). The injector shall be capable of accommodating at least three components in the sample loop, i.e. the injector programme shall be capable of controlling the “sandwich” injection of up to three components into the sample loop plus a solvent injection in one HPLC run.
- **Column:**
 - stainless steel;
 - inside diameter: 3 mm;
 - length: 250 mm;
 - temperature: 40 °C;
 - packing: reversed-phase C18 silica or equivalent;
 - particle size: 0,005 mm.

The resolution of the column shall be such that baseline separation of the components of interest is obtained.

The lifetime of the reversed-phase C18 column is very strongly influenced by the C18-silica bonding of the packing material. Therefore, columns equipped with a packing material containing monofunctional silanes with diisobutyl side-chain groups are preferred. These side groups sterically protect the key silanes from hydrolytic attack at low pH, making the stationary phase stable at such pH (pH 1).

- **UV detector:** wavelength 200 nm and 220 nm.

Additionally for determination of the linear oligomers:

- **Reagent pump**, including manometric module and pulse damper, giving a reagent flow rate of 0,25 ml/min and a reagent pressure drop of approximately 2 MPa.
- **Reaction coil:**
 - stainless steel or PEEK;
 - length: 3 m;

- inside diameter: 0,25 mm;
- temperature: 25 °C.
- **Fluorescence detector:**
 - excitation: 330 nm;
 - emission: 420 nm.

A schematic diagram of an HPLC apparatus is given in [Annex C](#).

6.2 Microbalance, accurate to 0,1 mg.

6.3 Ultrasonic vibration bath.

6.4 Volumetric flask, capacity 25 ml.

7 Test sample

For polyamide 6 and caprolactam, the maximum test sample size shall be 0,5 g. For mixtures of rearrangement products in water, a maximum test sample size of 1,25 g shall be used. The stability of the dissolved samples is limited, since caprolactam is hydrolysed to ϵ -aminocaproic acid. Do not keep for longer than 6 days.

8 Procedure

8.1 Calibration

8.1.1 Method A

To calibrate the column for the cyclic oligomers, prepare a series of at least three calibration solutions with concentrations increasing from 100 mg/l to 1 500 mg/l by dissolution in formic acid ([5.5](#)) for caprolactam ([5.6](#), $n = 1$), the cyclic dimer ([5.7](#), $n = 2$), the cyclic trimer ([5.8](#), $n = 3$) and the oligomer mixture ([5.9](#), $n = 4$ to $n = 6$). To calibrate the column for the linear oligomers, prepare a series of three calibration solutions with concentrations increasing from 2 mg/l to 15 mg/l by dissolving in formic acid ([5.5](#)) for ϵ -aminocaproic acid ([5.10](#)) and each of the linear oligomers from the dimer ([5.11](#)) to the hexamer ([5.15](#)).

If possible, adjust the concentrations of the calibration solutions in such way that the concentrations of the substance(s) in the sample solutions to be determined ([Clause 7](#)) are not in the lowest or highest part of the relevant calibration series.

Pump eluent A ([5.17](#)) through the column at a rate of 0,51 ml/min. Starting with the lowest concentration of a calibration series, inject 2 μ l of the calibration solution into the column in accordance with the injector programme given in [Annex A](#). Elute in accordance with the gradient timetable in [Annex A](#). Record the UV chromatogram. If applicable, immediately after UV detection, add the post-column reagent ([5.23](#)) at a rate of 0,25 ml/min, mixing the eluent and the reagent in the reaction coil. Record the fluorescence chromatogram. Measure the peak area of the component(s).

Repeat the calibration successively for the other calibration solutions of the same component and for the other calibration series.

NOTE Calibration can also be carried out using a commercially available PA6 polymer with a known concentration of cyclic and linear oligomers (see Reference [\[3\]](#)).

8.1.2 Method B

To calibrate the column for the caprolactam and its cyclic dimer, prepare a series of at least three calibration solutions with concentrations increasing from 100 mg/l to 1 500 mg/l by dissolution in formic acid (5.5) for caprolactam (5.6) and the cyclic dimer (5.7).

If possible, adjust the concentrations of the calibration solutions in such way that the concentrations of the substance(s) in the sample solution to be determined (Clause 7) are not in the lowest or highest part of the relevant calibration series.

Pump eluent B (5.17) through the column at a rate of 0,81 ml/min. Starting with the lowest concentration of a calibration series, inject 2 µl of the calibration solution into the column in accordance with the injector programme given in Annex B. Elute in accordance with the gradient timetable in Annex B. Record the UV chromatogram. Measure the peak area of the components.

Repeat the calibration successively for the other calibration solutions of the same component and for the other calibration series.

NOTE Calibration can also be carried out using a commercially available PA6 polymer with a known concentration of cyclic oligomers (see Reference [3]).


8.2 Determination

8.2.1 Method A

Introduce a test sample (see Clause 7) into a 25 ml volumetric flask. Add 20 ml of formic acid (5.5), close the flask and dissolve the sample, optionally by using ultrasonic vibration. Make up to the mark with formic acid (5.5). Pump eluent A (5.17) through the column at a flow rate of 0,51 ml/min. Inject 2 µl the sample solution into the column in accordance with the injector programme given in Annex A. Elute in accordance with the gradient timetable in Annex A. Record the UV chromatogram (see Figure 1). If applicable, immediately after UV detection, add the post-column reagent (5.23) at a flow rate of 0,25 ml/min, mixing the eluent and the reagent in the reaction coil. Record the fluorescence chromatogram (see Figure 2). Measure the peak areas of the cyclic oligomers (UV detection) and linear oligomers (fluorescence detection) up to and including the hexamer.


If necessary, e.g. if lower sensitivity is necessary due to high(er) concentrations, caprolactam can be detected by UV absorption at 220 nm instead of 200 nm.

NOTE The cyclic dimer of caprolactam elutes before caprolactam. See References [1], [2] and [3].



Key
 X retention time (min)
 Y response (mV)
 1 formic acid sandwich injection
 2 250 µl formic acid injection
 C1 to C6 cyclic oligomers
 L2 to L6 linear oligomers

Figure 1 — Determination of the cyclic oligomers of caprolactam — Method A, example — UV chromatogram



Key
 X retention time (min)
 Y response (mV)
 L1 to L6 linear oligomers


Figure 2 — Determination of the linear oligomers of caprolactam — Method A, example — Fluorescence chromatogram

8.2.2 Method B

Introduce a test sample (see [Clause 7](#)) into a 25 ml volumetric flask. Add 20 ml of formic acid ([5.5](#)), close the flask and dissolve the sample, optionally by using ultrasonic vibration. Make up to the mark with formic acid ([5.5](#)). Pump eluent B ([5.17](#)) through the column at a flow rate of 0,81 ml/min. Inject 2 µl of the sample solution into the column in accordance with the injector programme given in [Annex B](#). Elute in accordance with the gradient timetable in [Annex B](#). Record the UV chromatogram (see [Figure 3](#)). Measure the peak areas of the caprolactam and its dimer.

If necessary, e.g. if lower sensitivity is necessary due to high(er) concentrations, caprolactam can be detected by UV absorption at 220 nm instead of 200 nm.

NOTE The cyclic dimer of caprolactam elutes before caprolactam. See References [\[1\]](#), [\[2\]](#) and [\[3\]](#).



Key

X	retention time (min)
Y	response (mV)
1	formic acid sandwich injection
2	100 µl formic acid injection
C1 to C8	cyclic oligomers

Figure 3 — Determination of the cyclic oligomers of caprolactam — Method B, example — UV chromatogram

9 Calculations

9.1 Calculation of calibration factors

9.1.1 Calibration factors for cyclic oligomers

Calculate the calibration factor for the relevant component for each solution using [Formula \(1\)](#) (see Note 1):

$$f_c = \frac{A_c}{C_c} \quad (1)$$

where

- A_c is the peak area of the relevant component;
 C_c is the concentration of the relevant component in the calibration solution, in mg/l;
 f_c is the calibration factor for the relevant component.

For method A calculate the calibration factors for $n = 1$ up to and including $n = 3$ by taking the mean value of the three calculated calibration factors in the series of the relevant component. The calibration factors for the tetramer ($n = 4$) up to and including the hexamer ($n = 6$) are assumed to be constant. The calibration factor for the oligomers from $n = 4$ up to and including $n = 6$ is calculated from the chromatogram of the oligomer mixture, the calibration factors for $n = 1$ up to and including $n = 3$ being known. Calculate the calibration factor for $n = 4$ to $n = 6$ by taking the mean value of the three calculated calibration factors in the series.

For method B calculate the calibration factors by taking the mean value of the three calculated calibration factors.

The calibration shall be repeated if the relative difference between the calibration factors for two solutions is more than 5 %.

NOTE 1 The calibration factor for e.g. $n = 3$ can also be calculated by using for instance a calibration solution also containing caprolactam and/or the cyclic dimer if the latter are known.

NOTE 2 Usually, the linear oligomers are present in polyamide 6 at a concentration lower than that of the cyclic oligomers by approximately two orders of magnitude, and do not therefore interfere with the determination.

NOTE 3 The calibration factor can also be calculated by means of linear regression. The calibration curve can also be used for calculation of the relevant cyclic and/or linear oligomer concentrations.

9.1.2 Calibration factors for linear oligomers

Calculate the calibration factor for the linear oligomers for each solution using [Formula \(2\)](#):

$$f_{c1} = \frac{A_{c1}}{C_{c1}} \quad (2)$$

where

- A_{c1} is the peak area of the relevant component;
 C_{c1} is the concentration of the relevant component in the calibration solution, in mg/l;
 f_{c1} is the calibration factor for the relevant component.

Calculate the calibration factor for the linear components by taking the mean value of the three calculated calibration factors in the series. The calibration shall be repeated if the relative difference between the calibration factors for two solutions is more than 5 %.

NOTE The calibration factor can also be calculated by means of linear regression. The calibration curve can also be used for calculation of the relevant cyclic and/or linear oligomer concentrations.

9.2 Test sample

9.2.1 Calculation of cyclic oligomer concentrations

Calculate the concentration, in % by mass, of each of the cyclic oligomers from $n = 1$ up to and including $n = 6$ for method A or up to and including $n = 2$ for method B, using [Formula \(3\)](#) and rounding the result to 0,1 % for concentrations greater than or equal to 10 % and to 0,01 % for concentrations below 10 %:

$$C_{cs} = \frac{A_{cs} \times 2,5}{f_c \times m_g} \quad (3)$$

where

- A_{cs} is the peak area of the relevant component;
 f_c is the calibration factor for the relevant component, calculated in accordance with 9.1.1;
 m_g is the sample mass, in mg/25 ml;
 C_{cs} is the concentration of the relevant component in the sample, in % by mass.

NOTE For method A the concentration is calculated for each of the cyclic oligomers from $n = 1$ up to and including $n = 6$. For method B the concentration of caprolactam and its cyclic dimer is calculated.

9.2.2 Calculation of linear oligomer concentrations

Calculate the concentration, in mg/kg, of the linear oligomers from $n = 1$ up to and including $n = 6$ using Formula (4), rounding the result to 1 mg/kg in each case:

$$C_{cs} = \frac{A_{cs} \times 25\,000}{f_c \times m_g} \quad (4)$$

where

- A_{cs} is the peak area of the relevant component;
 f_c is the calibration factor for the relevant component, calculated in accordance with 9.1.2;
 m_g is the sample mass, in mg/25 ml;
 C_{cs} is the concentration of the relevant component in the sample, in mg/kg.

10 Precision

The reproducibility of this analytical method is not known because interlaboratory data are not available.

However, precision data within one laboratory have been determined by analysis of two reference polyamide 6 samples using method A and one polyamide 6 reference sample using method B, both over an extended period of time. The precision of this method is not expected to deviate significantly from the repeatability thus determined, the data for which is summarized in Tables 1 to 4 for method A and Table 5 for method B.

Table 1 — Repeatability data for determination of cyclic oligomers in a PA6 sample with a low cyclic oligomer content using method A

	$n = 1$	$n = 2$	$n = 3$	$n = 4$	$n = 5$	$n = 6$
Mean result, x (% by mass)	0,11	0,05	0,18	0,25	0,30	0,34
Standard deviation, s	0,007	0,008	0,008	0,01	0,01	0,02
Coefficient of variation (%)	6,9	15,6	4,5	3,9	4,4	6,9

Table 2 — Repeatability data for determination of cyclic oligomers in a PA6 sample with a high cyclic oligomer content using method A

	$n = 1$	$n = 2$	$n = 3$	$n = 4$	$n = 5$	$n = 6$
Mean result, x (% by mass)	8,37	0,33	0,50	0,47	0,40	0,34
Standard deviation, s	0,28	0,03	0,02	0,03	0,03	0,03
Coefficient of variation (%)	3,4	8,4	4,8	6,3	5,3	5,3

Table 3 — Repeatability data for determination of linear oligomers in a PA6 sample with a low linear oligomer content using method A

	$n = 1$	$n = 2$	$n = 3$	$n = 4$	$n = 5$	$n = 6$
Mean result, x (mg/kg)	<4	23	50	116	206	292

	<i>n</i> = 1	<i>n</i> = 2	<i>n</i> = 3	<i>n</i> = 4	<i>n</i> = 5	<i>n</i> = 6
Standard deviation, <i>s</i>	—	5	9	19	25	42
Coefficient of variation (%)	—	23	18	17	12	14

Table 4 — Repeatability data for determination of linear oligomers in a PA6 sample with a high linear oligomer content using method A

	<i>n</i> = 1	<i>n</i> = 2	<i>n</i> = 3	<i>n</i> = 4	<i>n</i> = 5	<i>n</i> = 6
Mean result, <i>x</i> (mg/kg)	141	59	82	148	203	252
Standard deviation, <i>s</i>	—	11	13	24	33	37
Coefficient of variation (%)	—	18	15	16	16	15

Table 5 — Repeatability data for determination of cyclic oligomers in a PA6 sample with a low cyclic oligomer content using method B

	<i>n</i> = 1	<i>n</i> = 2
Mean result, <i>x</i> (mg/kg)	0,13	0,03
Standard deviation, <i>s</i>	0,003	0,002
Coefficient of variation (%)	2,5	5,0

11 Test report

The test report shall include the following elements:

- a reference to this document, i.e. ISO 15033;
- all details necessary for complete identification of the sample, including type, manufacturer's code number, source, trade name, etc.;
- any deviation from the specifications for the HPLC equipment or from the procedure given in this document;
- any treatment of the sample prior to the analysis;
- the mass of the test sample used;
- the test method used A or B;
- the contents of the cyclic and linear oligomers of caprolactam, in % by mass and mg/kg, respectively;
- the relevant chromatograms.

Annex A (normative)

HPLC parameters and injector programme — Using eluent A

- a) Flow rate: 0,51 ml/min
- b) Solvent A: Eluent (5.17)
- c) Solvent B: Acetonitrile (5.4)
- d) Oven temperature: 40 °C
- e) Maximum pressure: 40 MPa
- f) Minimum pressure: Off
- g) Stop time: 47 min

Table A.1 — Eluent programme

Time min	Solvent A %	Solvent B %
0	99	1
36	68	32
36,1	99	1
47	99	1

Table A.2 — Injector programme

Line	Function	Volume µl	Time min	Solvent
1	Draw	3		Formic acid (5.5)
2	Draw	2		Sample solution
3	Draw	3		Formic acid (5.5)
4	Inject			
5	Wait		31	
6	Bypass			
7	Draw	250		Formic acid (5.5)
8	Inject			

NOTE Adjust slowdown of drawing, mixing and injection as appropriate and to suit the equipment.

Annex B (normative)

HPLC parameters and injector programme — Using eluent B

- a) Flow rate: 0,81 ml/min
- b) Solvent A: Eluent (5.17) method B
- c) Solvent B: Acetonitrile (5.4)
- d) Oven temperature: 40 °C
- e) Maximum pressure: 60 MPa
- f) Minimum pressure: Off
- g) Stop time: 15 min

Table B.1 — Eluent programme

Time min	Solvent A %	Solvent B %
0	90	10
5	68	32
5,1	90	10
12	90	10


Table B.2 — Injector programme

Line	Function	Volume µl	Time min	Solvent
1	Draw	3		Formic acid (5.5)
2	Draw	2		Sample solution
3	Draw	3		Formic acid (5.5)
4	Inject			
5	Wait		6	
6	Bypass			
7	Draw	100		Formic acid (5.5)
8	Inject			

NOTE Adjust slowdown of drawing, mixing and injection as appropriate and to suit the equipment.

Annex C (informative)

Schematic diagram of HPLC apparatus



Key

- | | | | |
|---|--------------|----|--|
| 1 | acetonitrile | 6 | UV detector (200 nm) |
| 2 | eluent | 7 | post-column reagent |
| 3 | eluent pump | 8 | reaction coil |
| 4 | injector | 9 | fluorescence detector (excitation 330 nm, emission 420 nm) |
| 5 | column | 10 | to waste |

Figure C.1 — Schematic diagram

Bibliography

- [1] MENERINK Y., PETERS R., KERKHOFF M., HELLENBRAND J., OMLOO H., ANDRIEN J. et al. Analysis of linear and cyclic oligomers in polyamide-6 without sample preparation by liquid chromatography using the sandwich injection method — I: Injection procedure and column stability. *J. Chromatogr. A.* 2000, **876** pp. 37–50
- [2] MENERINK Y., PETERS R., KERKHOFF M., HELLENBRAND J., OMLOO H., ANDRIEN J. et al. Analysis of linear and cyclic oligomers in polyamide-6 without sample preparation by liquid chromatography using the sandwich injection method — II: Methods of detection and quantification and overall long-term performance. *J. Chromatogr. A.* 2000, **878** pp. 45–55
- [3] MENERINK Y., PETERS R., VAN DER WAL S., CLAESSENS H.A., CRAMERS C.A. Analysis of linear and cyclic oligomers in polyamide-6 without sample preparation by liquid chromatography using the sandwich injection method — III: Separation mechanism and gradient optimization. *J. Chromatogr. A.* 2002, **949** pp. 307–326

British Standards Institution (BSI)

BSI is the national body responsible for preparing British Standards and other standards-related publications, information and services.

BSI is incorporated by Royal Charter. British Standards and other standardization products are published by BSI Standards Limited.

About us

We bring together business, industry, government, consumers, innovators and others to shape their combined experience and expertise into standards-based solutions.

The knowledge embodied in our standards has been carefully assembled in a dependable format and refined through our open consultation process. Organizations of all sizes and across all sectors choose standards to help them achieve their goals.

Information on standards

We can provide you with the knowledge that your organization needs to succeed. Find out more about British Standards by visiting our website at bsigroup.com/standards or contacting our Customer Services team or Knowledge Centre.

Buying standards

You can buy and download PDF versions of BSI publications, including British and adopted European and international standards, through our website at bsigroup.com/shop, where hard copies can also be purchased.

If you need international and foreign standards from other Standards Development Organizations, hard copies can be ordered from our Customer Services team.

Copyright in BSI publications

All the content in BSI publications, including British Standards, is the property of and copyrighted by BSI or some person or entity that owns copyright in the information used (such as the international standardization bodies) and has formally licensed such information to BSI for commercial publication and use.

Save for the provisions below, you may not transfer, share or disseminate any portion of the standard to any other person. You may not adapt, distribute, commercially exploit, or publicly display the standard or any portion thereof in any manner whatsoever without BSI's prior written consent.

Storing and using standards

Standards purchased in soft copy format:

- A British Standard purchased in soft copy format is licensed to a sole named user for personal or internal company use only.
- The standard may be stored on more than 1 device provided that it is accessible by the sole named user only and that only 1 copy is accessed at any one time.
- A single paper copy may be printed for personal or internal company use only.
- Standards purchased in hard copy format:
- A British Standard purchased in hard copy format is for personal or internal company use only.
- It may not be further reproduced – in any format – to create an additional copy. This includes scanning of the document.

If you need more than 1 copy of the document, or if you wish to share the document on an internal network, you can save money by choosing a subscription product (see 'Subscriptions').

Reproducing extracts

For permission to reproduce content from BSI publications contact the BSI Copyright & Licensing team.

Subscriptions

Our range of subscription services are designed to make using standards easier for you. For further information on our subscription products go to bsigroup.com/subscriptions.

With **British Standards Online (BSOL)** you'll have instant access to over 55,000 British and adopted European and international standards from your desktop. It's available 24/7 and is refreshed daily so you'll always be up to date.

You can keep in touch with standards developments and receive substantial discounts on the purchase price of standards, both in single copy and subscription format, by becoming a **BSI Subscribing Member**.

PLUS is an updating service exclusive to BSI Subscribing Members. You will automatically receive the latest hard copy of your standards when they're revised or replaced.

To find out more about becoming a BSI Subscribing Member and the benefits of membership, please visit bsigroup.com/shop.

With a **Multi-User Network Licence (MUNL)** you are able to host standards publications on your intranet. Licences can cover as few or as many users as you wish. With updates supplied as soon as they're available, you can be sure your documentation is current. For further information, email subscriptions@bsigroup.com.

Revisions

Our British Standards and other publications are updated by amendment or revision.

We continually improve the quality of our products and services to benefit your business. If you find an inaccuracy or ambiguity within a British Standard or other BSI publication please inform the Knowledge Centre.

Useful Contacts

Customer Services

Tel: +44 345 086 9001

Email (orders): orders@bsigroup.com

Email (enquiries): cservices@bsigroup.com

Subscriptions

Tel: +44 345 086 9001

Email: subscriptions@bsigroup.com

Knowledge Centre

Tel: +44 20 8996 7004

Email: knowledgecentre@bsigroup.com

Copyright & Licensing

Tel: +44 20 8996 7070

Email: copyright@bsigroup.com

BSI Group Headquarters

389 Chiswick High Road London W4 4AL UK